



CHM4B rabbit pAb

in-modifying protein 4b)) (SNF7-2) (hSnf7-2) (Vacuolar) (AA range: 50-100
) (SNF7-2) (hSnf7-2) (Vacuolar) (AA range: 50-100
) (SNF7-2) (hSnf7-2) (Vacuolar) (AA range: 50-100
) (SNF7-2) (hSnf7-2) (Vacuolar) (AA range: 50-100
) (SNF7-2) (hSnf7-2) (Vacuolar) (AA range: 50-100
) (SNF7-2) (hSnf7-2) (Vacuolar) (AA range: 50-100
n CHM4B
and 0.02% sodium azide.
m by affinity-chromatography
Peripheral membrane protein . clear envelope by CHMP7 zes transiently to the midbody 122861)
clear envelope by CHMP7 zes transiently to the midbody
clear envelope by CHMP7 zes transiently to the midbody 122861) eart and skeletal muscle. Also
2 y ay omornatogre

Nanjing BYabscience technology Co.,Ltd

网址: www.njbybio.com 官方热线: 025-5229-8998 监督电话: 15950492658







attention	Avoid repeated freezing and triawing:
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images